



Kanji Model® for suture training

JAPAN POLYMER TECHNOLOGY

Specification	
item#	Kanji Model WA
Size	W 90mm D 90mm H 3mm
Weight	35 g
resin	JPTC Polymer
Characters used	WA (和)
Oil for	non
Fixed type	Self-adhesive
Case for model	Available






Kanji Model® for suture training

JAPAN POLYMER TECHNOLOGY

Specification	
item#	Kanji Model NIPPON
Size	W 90mm D 90mm H 3mm
Weight	35 g
resin	JPTC Polymer
Characters used	NIPPON(日本)
Oil for	non
Fixed type	Self-adhesive
Case for model	Available



How to Use

JAPAN POLYMER TECHNOLOGY

Operating condition	<ul style="list-style-type: none"> • Flat place such as acrylic plate or glass surface • Rubber plates are not recommended as they generate gas and oil. • Surfaces with cork or powder will lose the stickiness of the model.
Confirmation before use	<ul style="list-style-type: none"> • Clean the surface on which the Kanji model is placed so that it is not contaminated with water, oil, dust, etc. • After using a model such as a 3D pad, wipe off the oil eluted from them before using the Kanji model.
Preparation for use	<ul style="list-style-type: none"> • Open the case for Kanji Model. image1 • Use the tip of the fingers to remove the edge of the model from the case. image2 • Support the peeled and sticky back side with a few fingertips. And softly peel off from the case and place in the dry box. Image3
Clean Up	<ul style="list-style-type: none"> • As with the preparation for use, slowly peel off the model with the tip of your finger and finish in the special case.



image1



image2



image3



Notes

JAPAN POLYMER TECHNOLOGY

Notes

Not for Medical Use

- If we hold the edge of the Kanji Model and try to peel it too strong, it will break.
- If we place the Kanji Model in a dirty place such as dust, it will lose its adhesive strength.
- If we apply a mold release agent, oil, acid, etc. to the Kanji Model, it will deteriorate.
- After using a model such as a 3D pad that elutes oil, be sure to wipe off the oil on the surface with alcohol before using.
- The rubber plate used in general dry boxes elutes gas and oil, which affects the Kanji model.
- If we place the Kanji Model on a document, the paper will not come off from it.
- The Kanji Model will deteriorate if placed in a hot car, such as in the summer.
- It will break if you pull it with more force than necessary in suture training.
- During suture training, pulling the thread with a strong force will break the Kanji Model.
- Be sure to clear off it in a case after use.
- After use, put the Kanji Model in a case and place it in a place out of direct sunlight.
- Please put the Kanji Model where the baby can't get at it



Our quality check items

JAPAN POLYMER TECHNOLOGY

Quality check items	•Check if the JPTC polymer is completely polymerized
	•Surface treatment is performed so that the Kanji Model surface does not become sticky.
	•Check the Kanji Model for defects and crevices.
	•Check the adhesive strength of the back of the Kanji Model on the front of the acrylic plate
	•Manufacture according to the JPTC manufacturing management control.
	•After manufacturing the Kanji Model, put it in a case and store it in the management box for 7 days.
	•Kanji Models in storage regularly open the case lid and inspect its quality.



Traidmark

JAPAN POLYMER TECHNOLOGY

Kanji Model is a trademark of JPTC

Traidmark : Kanji Model

Registration No.6304825

Registration Date : October 16, 2020

Kanji Model® is prohibited from being posted on the website or posted online without our permission.

Please find contact information here



JAPAN POLYMER TECHNOLOGY

Name : Masashi Inoue

E-Mail : inoue@jptc.co.jp