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Safety & Notes

Please read this instruction manual carefully and use the product safely and correctly.

Danger

Persons who use heart pacemaker or electro-medical apparatus may cause malpractice. Please do not use this product.

Warning

<p>When using this product, do not look directly into LED light.</p> <p>Do not use this product for any other purposes than making your face and body beautiful, or do not use the broken product. If you have any concerns about your skin and physical condition when using, please consult a doctor before using.</p> <p>For those who have any skin problems listed below, please do not use this product. It may worsen skin problems.</p> <ul style="list-style-type: none"> • Photosensitivity • Those who have skin disease or medical history of skin problems • Dermatitis (Atopic dermatitis, Contact dermatitis, excessive sunburn, etc.) • Metal allergy • Areas with purulence due to acne and pimple • Eczema, swelling, and warts (including virus) • Purpura • Melasma <p>Do not use this product on the skin and the areas that follow below. It may cause unexpected accidents or skin and body problems.</p> <ul style="list-style-type: none"> • Sensitive skin • Skin Paresthesia (especially those who cannot recognize heat) • Bugging mole • Cut • Surgical site due to accidents and orthopedic surgery • On implant (implant surgery metal / silicon / plastic, etc.) • Around the eyelids, eyeballs, and hollows around eyes • Adam's apple • Mouth and genital area, other mucous membrane • Thin skin area (immediately above the bone) • Around the heart <p>For those who have any disease listed below, please do not use this product. It may cause unexpected accidents or skin and body problems.</p> <ul style="list-style-type: none"> • Acute disease • Blood disease (hemorrhagic disease, hemophilia, etc.) • Those who are being under treatment or being followed up for abnormal blood pressure • Febrile disease • Malignant tumor • Epidemic (Infectious diseases) • Visceral disease • Facial teratopathy • Head and ophthalmologic sequelae due to traffic accidents • Those who have any abnormality in bone (osteoporosis, fracture, etc.) • Acute pain disease (sprain, muscle strain, etc.) <p>In the following cases, please stop using this product, turn the power off, and unplug the dedicated AC adapter from the outlet. There is a risk of fire or burns due to abnormal operation.</p> <ul style="list-style-type: none"> ● Smoke from this product ● Abnormal beating from this product ● Abnormal noise from this product ● Odor from this product <p>Never modify, disassemble or repair this product except by repair service engineers designated by us. There is a risk of fire or injury due to abnormal operation.</p> <p>Do not soak this product in water. There is a risk of electric shock or short circuit.</p> <p>Be sure to turn the power off when not in use or when cleaning this product. It may cause electric shock, short circuit, or fire.</p> <p>Unplug the dedicated AC adapter from the outlet except when charging. It may cause electric shock or earth leakage fire due to insulation deterioration.</p> <p>Do not connect or disconnect the dedicated AC adapter with wet hands. Also, make sure there is no water on the charging terminal before charging. There is a risk of electric shock or short circuit.</p> <p>This product contains a very strong infrared irradiation pulse program, infrared remote control devices such as television or air conditioner may not work. In the event that malfunction or any situation where remote control cannot be used arises, please turn off the power of this product.</p> <p>Keep this product out of reach of children. If not, it may cause accidents and problems.</p> <p>Do not use with pets or young children, and do not let them play with this product. If not, it may cause accidents or injury.</p> <p>For those who use this product for the first time or who have sensitive skin, shorten usage time to use it while checking the skin condition. If not, it may cause accidents and problems.</p> <p>Follow the instruction manual for usage time and frequency. If not, it may cause accidents or problems.</p> <p>If you feel pain or discomfort, stop using this product. If not, it may cause accidents and problems.</p> <p>Do not keep applying this product on the same area and be sure to move it on the skin. If not, it may cause burns or low temperature burns.</p>	<p>For those who have any of the following, please do not use this product. It may cause unexpected accidents or skin and body problems.</p> <ul style="list-style-type: none"> • Those who have epileptic seizures that cause to blinking lights • Allergic constitution • Those who have telangiectasia due to long-term use of steroid hormones or due to liver dysfunction • Those who are during menstruation • Those who are pregnant, one month after childbirth, and during nursing • Those who are under dental treatment, or temporomandibular disorders • Those who are in a growth process • Those who are after drinking alcohol or taking medicine (especially medicine that induces sleep) • Those who are very tired, who do not feel well, and who need rest • Those who are with body temperature of 38°C or higher (fever period) <p>Do not use this if there is any damage and abnormal heating in the cord and the dedicated AC adapter, or if it is loosely inserted into the outlet. There is a risk of fire, electric shock, short circuit, and burns.</p> <p>Do not hurt, damage, modify, forcibly bend, pull, twist, place a heavy object on the cord. There is a risk of fire, electric shock, short circuit, and burns.</p> <p>Be sure to use with AC 100-240V outlet. Also be sure to use the dedicated AC adapter for this product. If not, it may cause fire, electric shock, and failure.</p> <p>Remove precious metals such as earrings or necklaces for use. It may cause accidents or problems.</p> <p>Do not store this product in a bathroom or in a humid place, and do not place it on where is subject to water (such as a washbasin). There is a risk of fire or electric shock due to deterioration.</p> <p>Do not throw this product into fire or heat it, charge and use it at a high temperature, or leave it alone. It may cause heating and fire, burst, because this product has rechargeable battery built-in.</p> <p>Do not use this product near any flammable object (alcohol, thinner, benzine, liquid remover, spray, etc.) If not, it may cause fire or electric shock, or injury.</p>
--	---

Caution

<p>On the day of hair removal or shaving process, do not apply this product on the same area. There is a risk of swelling and redness.</p> <p>Do not drop, throw, and give a strong shock on this product. If not, it may cause accidents and problems.</p> <p>Do not place anything on this product.</p> <p>When unplugging the dedicated AC adapter, be sure to hold the adapter plug and not to pull on the cord.</p> <p>Not to store this product in the following places.</p> <ul style="list-style-type: none"> ● Place exposed to direct sunlight ● Place with high humidity ● Place with a lot of dust ● Place where there is a risk of ignition ● Place exposed to wind and rain ● Place at high temperature or below zero ● Place where chemicals or gas are generated ● Place with violent vibration <p>Do not wipe this product with chemicals such as alcohol, thinner, benzene, alkaline detergent, bleach, or disinfectant. If not, it may cause malfunction such as deformation, discoloration, crack, or damage.</p>	<p>For those who have any of the following, please do not use this product. It may cause unexpected accidents or skin and body problems.</p> <ul style="list-style-type: none"> • Those who have epileptic seizures that cause to blinking lights • Allergic constitution • Those who have telangiectasia due to long-term use of steroid hormones or due to liver dysfunction • Those who are during menstruation • Those who are pregnant, one month after childbirth, and during nursing • Those who are under dental treatment, or temporomandibular disorders • Those who are in a growth process • Those who are after drinking alcohol or taking medicine (especially medicine that induces sleep) • Those who are very tired, who do not feel well, and who need rest • Those who are with body temperature of 38°C or higher (fever period) <p>Do not use this if there is any damage and abnormal heating in the cord and the dedicated AC adapter, or if it is loosely inserted into the outlet. There is a risk of fire, electric shock, short circuit, and burns.</p> <p>Do not hurt, damage, modify, forcibly bend, pull, twist, place a heavy object on the cord. There is a risk of fire, electric shock, short circuit, and burns.</p> <p>Be sure to use with AC 100-240V outlet. Also be sure to use the dedicated AC adapter for this product. If not, it may cause fire, electric shock, and failure.</p> <p>Remove precious metals such as earrings or necklaces for use. It may cause accidents or problems.</p> <p>Do not store this product in a bathroom or in a humid place, and do not place it on where is subject to water (such as a washbasin). There is a risk of fire or electric shock due to deterioration.</p> <p>Do not throw this product into fire or heat it, charge and use it at a high temperature, or leave it alone. It may cause heating and fire, burst, because this product has rechargeable battery built-in.</p> <p>Do not use this product near any flammable object (alcohol, thinner, benzine, liquid remover, spray, etc.) If not, it may cause fire or electric shock, or injury.</p>
---	---

*Please keep this warranty in an accessible place for easy reference.
 *Misuse and incorrect storage of this product are not covered by the warranty, and repair is subject to charge even during the warranty period.

Before using

Frequency of use: Use within 24 minutes in total a day. Do not use esthetic function EMS/RF more than once a day.
 About skin: Before you use, be sure to perform a patch test (skin test) for each cosmetic you use. *This is a simple test to check any skin problems in advance.

- Patch test (skin test)
- 1 Wash the area (inside the arm) to be performed a patch test with soap, and cleanly wipe with towels.
 - 2 Apply cosmetics (tonic, beauty essence) inside part of the arm, and select the desincrustation (mode02) or iontophoresis (mode04).
 - 3 Apply the electrode of the product to the same area above, keep moving for about 3 minutes.
- If there is no abnormalities within 48 hours, you can use this product.
 * Be careful with any abnormalities on your skin when you use.
 * Do not use the product if it does not suit your skin condition.

How to operate

- Remove the main body from the charging stand.
- ① Power button/Start button
 Long press the button to turn the power on.
 Pressing the button briefly while the power is on will start each mode.
 * When you switch the mode while in use, select the mode again and press the start button briefly. Then new mode will start.
 * When you want to stop using, long press the button to turn off the power.
 - ② Mode button
 Press the button can select the mode (01→02→03→04→01). When you want to switch the mode while in use, select the mode again and press the power button/start button briefly to start new mode.
 - ③ Level button (EMS/RF)
 When you press the button when using mode 01 and 03, you can select strength of EMS/RF. You can also switch the strength during use.
 - ④ Function button
 When you press the button, you can select ON/OFF of esthetic function (EMS/desincrustation/RF/iontophoresis) in each mode. You can switch ON/OFF during use.
- *The end time may vary depending on the operating environment.

Tips for use

- How to hold the body
 Hold it the way that you can touch the ground electrode.
- Mode01 Texture · resilience
 - When you use EMS, make sure that the entire electrode part touches your skin, and slowly moving the body while using. (Reference : 2cm per second)
 - Mode02 Adolescent skin
 - Mode04 Bright skin
- When using desincrustation (mode 02) or iontophoresis (mode 04), use the cotton ring to attach the dedicated cotton.
- When using it with esthetic function OFF, the cotton is not necessary.
 After putting the cotton on, apply your own lotion or beauty essence to penetrate through entire cotton.
 After use, remove the cotton ring and dispose the used cotton, and wipe the LED lighting part and the electrode part of main body with a clean cloth. Wash the cotton ring with water and dry it well.
- Mode03 Rough skin · dry skin
- When using RF, use it that the electrode part touches the skin while moving the main body slowly. (Reference : 2cm per second)
- * Do not use around the eyelids or eyes.

After-sales service

- ◆ Rechargeable battery
 - Rechargeable battery is a consumable item. (Not covered by warranty)
 - If the operating time is extremely short after charging, it is considered that the rechargeable battery has reached the end of life.
 - If you want to replace the rechargeable battery, please contact the place of purchase or our customer service.
- ◆ How to dispose
 - When you dispose this product, please dispose it in accordance with your local regulations for sorting.
 - This product uses lithium ion rechargeable battery.

After-purchase servicing

- Warranty
 The warranty is attached to the back cover. Please fill in all necessary information such as location purchased and the date of purchase and keep it safe. Please read this warranty description carefully.
- Warranty period
 The guaranteed term is one year from the date of purchase.
- Request for repair
 Please contact our customer service if you still have problems after troubleshooting.
- Shipping costs shall be borne by the customer. We pay shipping costs after repair.
- During warranty period
 We will repair without charge according to the regulations of warranty. Please bring or send the product with this warranty and the proof-of-purchase, which can verify the date of purchase.
- After warranty period
 When the product still performs after repairing, we will repair with charge based on the needs of customers. Please contact the original distributor or our customer service.
- If any unsure about after service
 Please contact our customer service.

Please check this website for full text of the Manual / Warranty(English)