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Safety & Notes

Please read this instruction manual carefully and use the product safely and correctly.

Danger

Persons who use heart pacemaker or electro-medical apparatus may cause malpractice. Please do not use this product.

Warning

When using this product, do not look directly into LED light.

Do not use this product for any other purposes than making your face and body beautiful, or do not use the broken product. If you have any concerns about your skin and physical condition when using, please consult a doctor before using.

For those who have any skin problems listed below, please do not use this product. It may worsen skin problems.

- Those who have skin disease or medical history of skin problems
- Dermatitis (Atopic dermatitis, Contact dermatitis, excessive sunburn, etc.)
- Metal allergy
- Areas with purulence due to acne and pimple
- Eczema, swelling, and warts (including virus)
- Purpura
- Melasma

Do not use this product on the skin and the areas that follow below. It may cause unexpected accidents or skin and body problems.

- Sensitive skin

- Skin Paresthesia (especially those who cannot recognize heat)

- Bulging mole

- Cut

- Surgical site due to accidents and orthopedic surgery

- On implant/fimplant surgery metal / silicon / plastic, etc.)

- Around the eyelids, eyeballs, and hollows around eyes

- Adam's apple

- Mouth and genital area, other mucous membrane

- Thin skin area (immediately above the bone)

- Around the heart

For those who have any disease listed below, please do not use this product.

It may cause unexpected accidents or skin and body problems.

- Acute disease

- Blood disease (hemorrhagic disease, hemophilia, etc.)

- Those who are being under treatment or being followed up for blood pressure

- Fehile disease

- Malignant tumor

- Epidemic (Infection diseases)

- Visceral diseases

- Fever

- Head and ophthalmology sequelae due to traffic accidents

- Those who have any abnormality in bone (osteoporosis, fracture, etc.)

- Acute pain disease (sprain, muscle strain, etc.)

In the following cases, please stop using this product, turn the power off, and unplug the dedicated AC adapter from the outlet. There is a risk of fire or burns due to abnormal output.

●Smoke from this product ●Abnormal heating from this product ●Abnormal noise from this product ●Odor from this product

Never modify, disassemble or repair this product except by repair service engineers designated by us. There is a risk of fire or injury due to abnormal operation.

Do not soak this product in water. There is a risk of electric shock or short circuit.

Be sure to turn the power off when not in use or when cleaning this product. It may cause electric shock, short circuit, or fire.

Unplug the dedicated AC adapter from the outlet except when charging. It may cause electric shock or earth leakage fire due to insulation deterioration.

Do not connect or disconnect the dedicated AC adapter with wet hands. Also, make sure there is no water on the charging terminal before charging. There is a risk of electric shock or short circuit.

This product contains a very strong infrared irradiation pulse program, infrared remote control devices such as television or air conditioner may not work. In the event that malfunction or any situation where remote control cannot be used, please turn off the power of this product.

Keep this product out of reach of children. If not, it may cause accidents and problems.

Do not use with pets or young children, and do not let them play with this product. If not, it may cause accidents or injury.

For those who use this product for the first time or who have sensitive skin, shorten usage time to use it while checking the skin condition. If not, it may cause accidents and problems.

Follow the instruction manual for usage time and frequency. If not, it may cause accidents or problems.

If you feel pain or discomfort, stop using this product. If not, it may cause accidents and problems.

Do not keep applying this product on the same area and be sure to move it on the skin. If not, it may cause burns or low temperature burns.

Caution

On the day of hair removal or shaving process, do not apply this product on the same area. There is a risk of swelling and redness.

Do not drop, throw, and give a strong shock on this product. If not, it may cause accidents and problems.

Do not place anything on this product.

Not to store this product in the following places.

●Place exposed direct sunlight ●Face with high humidity ●Place with a lot of dust ●Place where there is a risk of ignition

●Place exposed to water or rain ●Face at high temperature or below zero ●Place where chemicals or gas are generated

●Place with violent vibration

Do not wipe this product with chemicals such as alcohol, thinner, benzene, alkaline detergent, bleach, or disinfectant. If not, it may cause malfunction such as deformation, discoloration, crack, or damage.

Before using

Frequency of use : Use within 20 minutes in total a day.

About skin : Before you use, be sure to perform a patch test (skin test) for each cosmetic you use.

This is a simple test to check any skin problems in advance.

- Patch test (skin test)**
- 1 Wash the area (inside the arm) to be performed a patch test with soap, and cleanly wipe with towels.
 - 2 Apply cosmetics (lotion, beauty essence) inside part of the arm, and select the desincrustation(P05) or ioncleanse(P06).
 - 3 Apply the Ion head to the same area above, keep moving for about 3 minutes.
- If there is no abnormalities within 48 hours, you can use this product.
- * Be careful with any abnormalities on your skin when you use.
- * Do not use the product if it does not suit your skin condition.

How to use

Program / Remaining battery indication on display

| | | |
|--|------|---|
| | 100% | In standby mode, the value of remaining battery blinks. |
| | 75% | Lights up while each LED program (P01 to P04) is running. |
| | 50% | Blinks one by one while ion program (P05 / P06) is running. |
| | 25% | When turn off the power, the value of remaining battery indicates for 1 second. |

Remaining battery display

①Put your finger (index finger) through the hold ring to hold the main body so that you can touch the hand electrode with hand.
②Turn the power on and select the program.

-Long press the power button to power on and enter standby mode. (startup sound → standby mode:remaining battery indicated)
-By a short press, you can switch LED program (P01 to P04) or ion program (P05 / P06).

-When you press the operation button while each program is running, vibration function (interlocking mode) will activate.
※ Long press the power button to power off both in standby mode and during each program operation. (power off sound → remaining battery indicated for 1 second → power off)

※ Auto-off function built-in: When the program ends, the power turns off automatically.

In addition, if there is no operation for 60 seconds in standby mode, the power turns off automatically. (power off sound → remaining battery indicates for 1 second → power off)

●LED program (P01 to P04)

Please select the program according to your skin condition and troubles.
For the details of each program, please refer to the program list .

●Desincrustation program (P05)

Please use it with your skin wet.
After the treatment, please wipe off with cotton cloth soaked with lotion, or rinse off lightly.

●Ionophoresis program (P06)

Please use it with lotion / beauty essence (contains vitamin C, vitamin C derivative).
-Put the ion head and LED lighting part on the skin and move it slowly.

-Touch sensor is equipped, please keep the ion head touching the skin while using it.

-When the program (about 5 minutes) ends, the power-off sound is emitted and the power turns off automatically.

【Ion program】

-Press the button in standby mode or during each LED program operation switches to ion program, and desincrustation (P05) starts automatically.

-Each time you press the button, ion program (P05/P06) can be switched.

【Vibration function】

-Pressing the button while each program is running activates vibration function (interlocking mode).

-Each time you press the button, you can switch the Interlocking mode → weak → strong → OFF.

※ If the program is changed, vibration function will be automatically turned off. please set it again.

※ You cannot use vibration function without program.

※ P02 will not activate the vibration function because they cannot be used together in the program.

④Care for the main body after using and store it in the dedicated porch.

After-sales service

Rechargeable battery

◆Rechargeable battery
-Rechargeable battery is a consumable item. (Not covered by warranty)

-If the operating time is extremely short after charging, it is considered that the rechargeable battery has reached the end of life.
-If you want to replace the rechargeable battery, please contact the place of purchase or our customer service.

How to dispose

-When you dispose this product, please dispose it in accordance with your local regulations for sorting.

-This product uses lithium ion rechargeable battery.

After-purchase service

Warranty

The warranty is attached to the back cover. Please fill in all necessary information such as location purchased and the date of purchase and keep it safe.

Please read this warranty description carefully.

Warranty period

The guaranteed term is one year from the date of purchase.

Request for repair

Please contact our customer service if you still have problems after troubleshooting.

※ Shipping costs shall be borne by the customer. We pay shipping costs after repair.

During warranty period

We will repair without charge according to the regulations of warranty. Please bring or send the product with this warranty and the proof-of-purchase, which can verify the date of purchase.

After warranty period

We will repair with charge based on the needs of customers. Please contact the original distributor or our customer service.

If any unsure about after service

Please contact our customer service.

Please check this website for full text of the Manual / Warranty(English).